



Da Vinci® Instruments

## **Force Bipolar**

With DualGrip  
technology

Added control.\*  
Nonstop progress.\*\*

\* Compared to other da Vinci  
8 mm instruments.  
Data on file.

\*\* Nonstop progress refers  
to reduced instrument  
exchanges.

DAVINCI™ by Intuitive





## Added control.\* Nonstop progress.\*\*

Get a new level of grasping control throughout your case\* with Force Bipolar and DualGrip technology.

Close surgeon collaboration aimed at increasing efficiency in da Vinci® cases led to a design that features two distinct grip modes.

Stay focused on the surgical field and shift grip strength when you need to, instead of asking your assistant to switch in another grasper.

With the Force Bipolar in your toolkit, you can streamline setup and minimize workflow disruptions by reducing instrument exchanges.

\* Compared to other da Vinci 8 mm instruments.  
Data on file. Available for the da Vinci X® and/or da Vinci Xi® surgical systems

\*\* Nonstop progress refers to reduced instrument exchanges.

# Take grasping control to a new level.\*

Increase the strength of your grip seamlessly with DualGrip technology.

Simply switch the mode from "Default" to "Strong" using your foot pedal to get a stronger hold of the suture or tissue you need to handle, instead of interrupting the procedure to bring in a new instrument.



\* Compared to other da Vinci 8 mm instruments.

# Be efficient.

Force Bipolar combines the abilities of two grasping instruments, which enables you to increase efficiencies<sup>1</sup> and your hospital to potentially save on instrument costs.<sup>2</sup>


## Now you can:

Simplify instrument selection with a grasper that covers multiple needs.

Minimize disruptions to your workflow by reducing instrument exchanges.

Potentially save time during your cases, delivering possible cost savings to your hospital.<sup>2</sup>

## Potential savings from reduction in number of instruments.<sup>2</sup>



1. Due to DualGrip technology and minimized instrument exchanges

2. Internal analysis on file.



### Da Vinci® Force Bipolar

1 unit (10 uses) | 470405

da Vinci X and da Vinci Xi

#### You can take advantage of Force Bipolar for:

grasping

retraction

dissection

bipolar coagulation of tissue

#### To contact a representative or receive additional information:

Visit [www.intuitive.com](http://www.intuitive.com) or call Intuitive Surgical Customer Service:  
in Europe at +41.21.821.2020 or +800.0821.2020 (toll free)

The EndoWrist® Force Bipolar is intended to be used with compatible systems for endoscopic manipulation of tissue, including dissection, grasping, retraction, and bipolar coagulation of tissue.

The EndoWrist® Force Bipolar is a class IIb medical device CE marked (CE0543) under the European Medical Devices Directive (93/42/EEC) manufactured by Intuitive Surgical, Inc. Refer to Instructions For Use before use.



DAVINCI™ by Intuitive